



# Recombinant Human PCNA-interacting partner (PARPBP)

<b>Product Code</b>	CSB-EP868339HU
<b>Abbreviation</b>	PARPBP
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9NWS1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MAVFNQKSVS DMIKEFRKNW RALCNSERTT LCGADSMMLLA LQLSMAENNK QHSGEFTVSL SDVLLTWKYL LHEKLNLPVE NMDVTDHYED VRKIYDDFLK NSNMLDLIDV YQKCRALTSN CENYNTVSPS QLLDFLSGKQ YAVGDETDLS IPTSPTSKYN RDNEKVQLLA RKIIFSYNL LVNSKNDLAV AYILNIPDRG LGREAFTDLK HAAREKQMSI FLVATSFIRT IELGGKGYAP PPSDPLRTHV KGLSNFINFI DKLDEILGEI PNPSIAGGQI LSVIKMQLIK GQNSRDPFCK AIEEVAQDLD LRIKNIINSQ EGVVALSTTD ISPARPKSHA INHGTAFCGR DTVKALLVLL DEEAANAPTK NKAELLYDEE NTIHHHGTSI LTLFRSPTQV NNSIKPLRER ICVSMQEKKI KMKQTLIRSQ FACTYKDDYM ISKDNWNNVN LASKPLCVLY MENDLSEGVN PSVGRSTIGT SFGNVHLDRS KNEKVSRSKST SQTGNKSSKR KQVDLDGENI LCDNRNEPPQ HKNKIPKKS NDSQNRLYGK LAKVAKSNKC TAKDKLISGQ AKLTQFFRL
<b>Source</b>	E.coli
<b>Target Names</b>	PARPBP
<b>Protein Names</b>	Recommended name: PCNA-interacting partner Short name= PARI Alternative name(s): PARP-1 binding protein PARP1-binding protein Short name= PARPBP
<b>Expression Region</b>	1-579
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.