



# Recombinant Human Serine/threonine-protein kinase 32B (STK32B)

<b>Product Code</b>	CSB-YP868353HU
<b>Abbreviation</b>	STK32B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9NY57
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MGGNHSKPP VFDENEEVNF DHFQILRAIG KGSFGKVCIV QKRDTKKMYA MKYMNKQKCI ERDEVNRVFR ELQIMQGLEH PFLVNLWYSF QDEEDMFMVV DLLLGGDLRY HLQQNVHFTG GTVKLYICEL ALALEYLQRY HIIHRDIKPD NILLDEHGHV HITDFNIATV VKGAERASSM AGTKPYMAPE VFQVYMDRGP GYSYPVDWWS LGITAYELLR GWRPYEIHVS TPIDEILNMF KVERVHYSST WCKGMVALLR KLLTKDPESR VSSLHDIQSV PYLADMNWDVA VFKKALMPGF VPNKGRLNCD PTFELEEMIL ESKPLHKKKK RLAKNRSRDG TKDSCPLNGH LQHCLETVRE EFIIFNREKL RRQQGQGSQ LDTDSRGGGQ AQSKLQDGCN NLLTHTCTR GCSS
<b>Source</b>	Yeast
<b>Target Names</b>	STK32B
<b>Protein Names</b>	Recommended name: Serine/threonine-protein kinase 32B EC= 2.7.11.1 Alternative name(s): YANK2
<b>Expression Region</b>	1-414
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.