



Recombinant Human U3 small nucleolar RNA-associated protein 6 homolog (UTP6)

Product Code	CSB-MP868357HU
Abbreviation	UTP6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NYH9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	AEIIQERIE DRLPELEQLE RIGLFSHAEI KAIKKASDL EYKIQRRTLF KEDFINYVQY EINLLELIQR RRTRIGYSFK KDEIENSIVH RVQGVFQRAS AKWKDDVQLW LSYVAFCKKW ATKTRLSKVF SAMLAIHSNK PALWIMAAKW EMEDRLSSES ARQLFLRALR FHPECPKLYK EYFRMELMHA EKLKKEKEEF EKASMDVENP DYSEEILKGE LAWIIYKNSV SIIKGAEFHV SLLSIAQLFD FAKDLQKEIY DDLQALHTDD PLTWDYVARR ELEIESQTEE QPTTKQAKAV EVGRKEERCC AVYEEAVKTL PTEAMWKCYI TFCLERFTKK SNSGFLRGKR LERTMTVFRK AHELKLLSEC QYKQLSVSLL CYNFLREALE VAVAGTELFR DSGTMWQLKL QVLIKESKSPD IAMLFEAFV HLPKQVCLPL WISWAEWSEG AKSQEDTEAV FKKALLAVIG ADSVTLKNKY LDWAYRSGGY K KARAVFKSL QESRPFSVDF FRKMIQFEKE QESCNMANIR EYYERALREF GSADSDLWMD YMKEELNHPL GRPENCGQIY WRAMKMLQGE SAEAFVAKHA MHQTGHL
Source	Mammalian cell
Target Names	UTP6
Protein Names	Recommended name: U3 small nucleolar RNA-associated protein 6 homolog Alternative name(s): Hepatocellular carcinoma-associated antigen 66 Multiple hat domains protein
Expression Region	2-597
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.