



Recombinant Human TGF-beta-activated kinase 1 and MAP3K7-binding protein 2 (TAB2), partial

Product Code	CSB-MP868358HU
Abbreviation	TAB2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NYJ8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MRNQPTLFISTNSGASAASRNMSGQVSMGPAFIHHHPPKSRAIGNNSATSPR VVVTQPNTKYTFKITVSPNKPPAVSPGVVSPTELTNLLNHPDHYVETENIQHL TDPTLAHVDRISETRKLMSGDDAAYTQALLVHQKARMERLQRELEIQKKKLD KLKSEVNEMENNLTRRRLKRSNSISQIPSLEEMQQLRSCNRQLQIDIDCLTKEI DLFQARGPHFNPSAIHNFYDNIGFVGPVPPKPKDQRSIIKTPKTQDTEDEGAQ WNCTACTFLNHPALIRCEQCEMPRHF
Source	Mammalian cell
Target Names	TAB2
Protein Names	Recommended name: TGF-beta-activated kinase 1 and MAP3K7-binding protein 2 Alternative name(s): Mitogen-activated protein kinase kinase kinase 7-interacting protein 2 TAK1-binding protein 2 Short name= TAB-2 TGF-beta-activated kinase 1-bi
Expression Region	401-693
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.