



Recombinant Human Pleckstrin-2 (PLEK2)

Product Code	CSB-YP868361HU
Abbreviation	PLEK2
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q9NYT0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MEDGVLKEGF LVKRGHIVHN WKARWFILRQ NTLVYYKLEG GRRVTPPKGR ILLDGCTITC PCLEYENRPL LIKLTQTST EYFLEACSRE ERDAWAFEIT GAIHAGQPGK VQQLHSLRNS FKLPPHISLH RIVDKMHDSN TGIRSSPNME QGSTYKKTFL GSSLVDWLIS NSFTASRLEA VTLASMLMEE NFLRPVGVRS MGAIRSGDLA EQFLDDSTAL YTFAESYKKK ISPKEEISLS TVELSGTVVK QGYLAKQGHK RKNWKVRRFV LRKDPAFLHY YDPSKEENRP VGGFSLRGS LSALEDNQVP TGVKGNVQGN LFKVITKDDT HYYIQASSKA ERAEWIEAIK KLT
Source	Yeast
Target Names	PLEK2
Protein Names	Recommended name: Pleckstrin-2
Expression Region	1-353
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>