



CXXC-type zinc finger protein 1 (CXXC1)

Product Code	CSB-EP868380HU-B
Abbreviation	CXXC1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9P0U4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MEGDGSDPEP PDAGEDSKSE NGENAPIYCI CRKPDINCFM IGCNDCNEWF HGDCIRITEK MAKAIREWYC RECREKDPKL EIRYRHKKS ERDGNERDSS EPRDEGGGRK RPVPDPDLQR RAGSGTG VGA MLARGSASPH KSSPQPLVAT PSQHHQQQQQ QIKRSARMCG ECEACRRTED CGHCDFCRDM KKFGGPNKIR QKCRLRQCQL RARESYKYFP SSLSPVTPSE SLPRRRPLP TQQQPQPSQK LGRIREDEGA VASSTVKEPP EATATPEPLS DEDLPLDPDL YQDFCAGAFD DHGLPWMSDT EESPFLDPAL RKRAVKVKHV KRREKKSEKK KEERYKRHRQ KQKHKDKWKH PERADAKDPA SLPQCLGPGC VRPAQPSSKY CSDDCGMKLA ANRIYEILPQ RIQQWQQSPC IAEEHGKLL ERIRREQQSA RTRLQEMERR FHELEAILR AKQQAVREDE ESNEGDSDDT DLQIFCVSCG HPINPRVALR HMERCYAKYE SQTSGSMYP TRIEGATRLF CDVYNPQSKT YCKRLQVLCPEHSRDPKVPA DEVCGCPLVR DVFELTGDFC RLPKRQCNRH YCWEKLRAE VDLERVRVWY KLDELFEQER NVRTAMTNRA GLLALMLHQT IQHDPLTTDL RSSADR
Source	E.coli
Target Names	CXXC1
Protein Names	Recommended name: CpG-binding protein Alternative name(s): CXXC-type zinc finger protein 1 PHD finger and CXXC domain-containing protein 1
Expression Region	1-656
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself.
Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.