



Recombinant Human Putative serine/threonine-protein phosphatase 4 regulatory subunit 1-like (PPP4R1L)

Product Code	CSB-MP868382HU
Abbreviation	PPP4R1L
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9P1A2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAEIPLYFVD LQDDLDDYGF EDYGTDCDNM RVTAFLDIPG QDNLPLTRL EKYAFSENTF NRQIIARGLL DIFRDFGNNE EDFLTVM EIV VRLSEDAEPT VRTELMEQIP PIAIFLQENR SNFPVVLSEY LIPIVVRYLT DPNNQIICKM ASMLSKSTVE RLLLPRFC EL CGDRKLFQVR KVCAANFGDI CHAVGQEATE KFLIPKFFEL CSDAVWGM RK ACAECFTAVS HSSSPGVRRT QLFPLFIRLV SDPCRWWHQA AFQSLGPFTS TFANPSRAGL YLREDGALSI WPLTQDLDSG FASGSPAPSS GGNTSPASLT RSKAPVRSEP ELPVEGTS AK TSDCPHSSSS SDGPAESPVE SCVSAGAEWT RVSPETSACS KLSMDNDLPI SSYPGSDSWA CPGNTEDVFS HFLYW
Source	Mammalian cell
Target Names	PPP4R1L
Protein Names	Recommended name: Putative serine/threonine-protein phosphatase 4 regulatory subunit 1-like
Expression Region	1-415
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.