



# Recombinant Human Succinate--CoA ligase [ADP-forming] subunit beta, mitochondrial (SUCLA2)

<b>Product Code</b>	CSB-MP868398HU
<b>Abbreviation</b>	SUCLA2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9P2R7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	LSLHEYMS MELLQEAGVS VPKGYYAKSP DEAYAIKKL GSKDVIKAAQ VLAGGRGKGT FESGLKGGVK IVFSPEEAKA VSSQMIGKKL FTKQTGEKGR ICNQVLVCER KYPRREYYFA ITMERSFQGP VLGSSHGGV NIEDVAAESP EAIKEPIDI EEGIKKEQAL QLAQKMGFPP NIVESAAENM VKLYSLFLKY DATMIEINPM VEDSDGAVLC MDAKINFDN SAYRQKKIFD LQDWTQEDER DKDAAKANLN YIGLDGNIGC LVNGAGLAMA TMDIILHGG TPANFLDVGG GATVHQVTEA FKLITSDKKV LAILVNIFGG IMRCDVIAQG IVMVAVKLEI KIPVVVRLQG TRVDDAKALI ADSGLKILAC DDLDEAARMV VKLSEIVTLA KQAHVDVKFQ LPI
<b>Source</b>	Mammalian cell
<b>Target Names</b>	SUCLA2
<b>Protein Names</b>	Recommended name: Succinyl-CoA ligase [ADP-forming] subunit beta, mitochondrial EC= 6.2.1.5 Alternative name(s): ATP-specific succinyl-CoA synthetase subunit beta Renal carcinoma antigen NY-REN-39 Succinyl-CoA synthetase beta-A chain S
<b>Expression Region</b>	53-463
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.