



Recombinant Schizosaccharomyces pombe Putative Rho GTPase-activating protein C1565.02c (SPAC1565.02c)

Product Code	CSB-EP868404SXV-B
Abbreviation	SPAC1565.02c
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9P3B1
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MEHDLERLCF IGGYDNDNDK VIVVVTKNLE LFKKYDDINL IKEAYNHVHK LIQKDERYTA VFFAHDSTVF SYLGLSLKAY YGMDYYLHKN VKAVYVIHTD WMSKVAIRTL LSIASPKFTR KFRYLNSISD LNKYIPLSHL KLPPIVYEFD RNVEPAIFPS NPSATSNFIF PVQQWTRPPD ALVEGMQVVS KSLQTEGLFR KSCSRKHLDI VIELYDNGCM VDLEAFGPIA ACLIKHLFR SLPSPLFSAE FLNGLTDHMD SGIDYAVSLQ KLIDASMDKN SQKLARLIFS LLYQITQHEQ ENMMNAQNLA LCIGPSFSKA DNIAELMSMK DKEYNPYCYF LEYAILHWNT LFANEENWDS YIPQDPTLLL PKTP
Source	E.coli
Target Names	SPAC1565.02c
Protein Names	Recommended name: Putative Rho GTPase-activating protein C1565.02c
Expression Region	1-374
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.