



Recombinant *Schizosaccharomyces pombe* Uncharacterized protein C21C3.11 (SPBC21C3.11)

Product Code	CSB-MP868449SXV
Abbreviation	SPBC21C3.11
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9P7L2
Product Type	Recombinant Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MATIYACRGF HRVPVKLSPS STLQEVLSS YKQLGFSDWH NLELLHGDKK VDTSLLLRLS GIINGAKLIV KESATNQSGK SSSSISPQSK KIKVALQLPG AARIIDEASS ETSIKQLLER HSLLTkVSHV LINGRNFKSE EFDNPLLLYG IREGSILIRL FPIKAQQSIV SEQAPVSQTF NGDVKEKKNA DLMDIESENK KDDIVESFPK YPVD AHLLE PLPTPIPSLP STPSSYQNL P SQSLTGESLP TVSNQEKDEG VIEKVAVNNT PSVSSKSPFP KKKSFSSMLA QVKKEKAENN GSDGYDLQPT KSQLELYQSI LRKRANQVSS TSLTKSSSPK PLPSSAIVKF DFGNGKSIVH EFSKDDNIET LRAVASHLS PEESTSFQLT FSNYEALPTT GLIVEHIGRA VVRVHTISDP VYAQR
Source	Mammalian cell
Target Names	ubx4
Protein Names	Recommended name: Uncharacterized protein C21C3.11
Expression Region	1-425
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.