



# Recombinant Candida albicans Lipase 5 (LIP5)

<b>Product Code</b>	CSB-BP868468CZD
<b>Abbreviation</b>	LIP5
<b>Storage</b>	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
<b>Uniprot No.</b>	Q9P8W0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Candida albicans (strain SC5314 / ATCC MYA-2876) (Yeast)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	GLIFPT KPSSDPFYNP PKGFENAAVG DILQSRATPK SITGGFTPLK IQNSWQLLVR SEDSFGNPNV IVTTVIEPVN ADPSKIASYQ VFEDAAKADC APSYALQFGS DLTTFVTQAE MYLMAPLLDQ GYYVVSVDYE GPKSTFTIGK QSGQAVLNSI RAALKSGKIT NIKDDAKVVM WGYSGGSLAS GWAAALQPSY APELGGNLLG AALGGFVTNI TATAQATDGT VFAGIVANAL GGVANEYPEF KSILQSDTDK KSVFDEFDGH CLIDGVLNYI GTSFLTGDHK IFKTGWDILK NPKIGKVVED NGLVYQKQLV PKIPVFIYHG SIDQIVPIVD TKKTYQNWCD AGISSLEFAE DASNGHLTEA IMGAPAALTW IIDRFDGKQT VSGCQHIQRF SNLEYPNIPS SIANYFKAAM DVVLHLGLGP DVQKDQVSPE GIKKLGSIEM RWL
<b>Source</b>	Baculovirus
<b>Target Names</b>	LIP5
<b>Protein Names</b>	Recommended name: Lipase 5 EC= 3.1.1.3
<b>Expression Region</b>	15-463
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>