



# Recombinant Human Insulin-like growth factor-binding protein complex acid labile subunit (IGFALS)

<b>Product Code</b>	CSB-MP011094HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P35858
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	ADP GTPGEAEGPA CPAACVCSYD DDADELSVFC SSRNLTRLPD GVPGGTQALW LDGNNLSSVP PAAFQNLSSL GFLNLQGGQL GSLEPQALLG LENLCHLHLE RNQLRSLALG TFAHTPALAS LGLSNNRLSR LEDGLFEGLG SLWDLNLGWN SLAVLPDAAF RGLGSLRELV LAGNRLAYLQ PALFSGLAEL RELDLSRNAL RAIKANVFVQ LPRLQKLYLD RNLIAAVAPG AFLGLKALRW LDLSHNRVAG LLEDTFPGLL GLRVLRLSHN AIASLRPRTF KDLHFLEELQ LGHNRIKRLA ERSFEGLGQL EVLTLDHNLQ QEVKAGAFGLG LTNVAVMNLN GNCLRNLPDQ VFRGLGKLHS LHLEGSCLGR IRPHTFTGLS GLRRLFLKDN GLVGIEEQSL WGLAELELD LTSNQLTHLP HRLFQGLGKL EYLLSRNRL AELPADALGP LQRAFWLDVS HNRLEALPNS LLAPLGLRLRY LSLRNNSLRT FTPQPPGLER LWLEGNPWDC GCPLKALRDF ALQNPSAVPR FVQAICEGDD CQPPAYTYNN ITCASPPEVV GLDLRDLSEA HFAPC
<b>Source</b>	Mammalian cell
<b>Target Names</b>	IGFALS
<b>Protein Names</b>	Recommended name: Insulin-like growth factor-binding protein complex acid labile subunit Short name= ALS
<b>Expression Region</b>	28-605
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This protein is a serum protein that binds insulin-like growth factors, increasing their half-life and their vascular localization. Production of the encoded protein, which contains twenty leucine-rich repeats, is stimulated by growth hormone. Defects in this gene are a cause of acid-labile subunit deficiency, which manifests itself in a delayed and slow puberty. Three transcript variants encoding two different isoforms have been found for this gene.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.