



# Recombinant Chicken Estrogen receptor beta (ESR2)

<b>Product Code</b>	CSB-BP868699CH
<b>Abbreviation</b>	ESR2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9PTU5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Gallus gallus (Chicken)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MAFCSPAMLN YNIASNFGDS ESASVRQTSS PSQLWSAPGH LSPLTLHCQS SLLYAEQPKS PWCEVRPLDP VLPVTRETLK RKTNGSDCTS PIASNPGSKR DAHFCVCS D YASGYHYGVW SCEGCKAFFK RSIQGHNDYI CPATNQCTID KNRRKSCQAC RLRKCYEVGM MKCGSRRERC GYRILRRHRN SEDCMGKTKK YNEAATRVKE ILLSTVSPEQ FVLTLLEAEP PNVLVS RPSK PFTEASMMMS LTKLADKELV HMIGWAKKIP GFIDLSLYDQ VRLLESCWME VLMIGLMWRS IDHPGKLIFA PDLVLDRDEG KCVEGILEIF DMLLAMTSRF RELKLQHKEY LCVKAMILLN SSMFPLSPEE PESNRKLHHL LNVVTDALVW VIAKSGIPSQ QQTTRLANLL MLLSHVRHAS NKGMEHLLSM KCKNVVPVYD LLEMLNAHT LRGQRKSPVT HPEFEQVSHF QV
<b>Source</b>	Baculovirus
<b>Target Names</b>	ESR2
<b>Protein Names</b>	Recommended name: Estrogen receptor beta Short name= ER-beta Short name= cERb Alternative name(s): Nuclear receptor subfamily 3 group A member 2
<b>Expression Region</b>	1-472
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	This gene encodes a member of the family of estrogen receptors and superfamily of nuclear receptor transcription factors. The gene product contains an N-terminal DNA binding domain and C-terminal ligand binding domain and is localized to the nucleus, cytoplasm, and mitochondria. Upon binding to 17beta-estradiol or related ligands, the encoded protein forms homo- or hetero-dimers that interact with specific DNA sequences to activate transcription. Some isoforms dominantly inhibit the activity of other estrogen receptor family



members. Several alternatively spliced transcript variants of this gene have been described, but the full-length nature of some of these variants has not been fully characterized.

---

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.