



Recombinant Salmonella typhimurium Secreted effector protein SseD (sseD)

Product Code	CSB-MP870836SXB
Abbreviation	sseD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9R803
Product Type	Recombinant Protein
Immunogen Species	Salmonella typhimurium (strain LT2 / SGSC1412 / ATCC 700720)
Purity	>85% (SDS-PAGE)
Sequence	MEASNVALVL PAPSLLTPSS TPSPSGEGMG TESMLLLFDD IWMKLMELAK KLRDIMRSYN VEKQRLAWEL QVNVLTQTMK TIDEAFRASM ITAGGAMLSG VLTIGLGAVG GETGLIAGQA VGHTAGGVMG LGAGVAQRQS DQDKAIADLQ QNGAQSYNKS LTEIMEKATE IMQQIIGVGS SLVTVLAEIL RALTR
Source	Mammalian cell
Target Names	sseD
Protein Names	Recommended name: Secreted effector protein SseD Alternative name(s): Secretion system effector D
Expression Region	1-195
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.