



Recombinant Arabidopsis thaliana Probable ADP-ribosylation factor At2g15310 (At2g15310)

Product Code	CSB-YP871067DOA
Abbreviation	At2g15310
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9SHU5
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	GARFSRIAK RFLPKSKVRI LMVGLDGSGK TTILYKCLKG EVVTTVPTIG FNLETVEYKG INFTVWDIGG QEKIRKLWRH YFQNAQGLIF VVDSSDSERL SEARNELHRI LTDNELEGAC VLVFANKQDS RNALPVAEVA NKLGLHLSLK RCWLIQGTSA ISGQGLYEGL EWLSTTIPNK PERSTSVSSF RSDSYERKLV RGPRY
Source	Yeast
Target Names	At2g15310
Protein Names	Recommended name: Probable ADP-ribosylation factor At2g15310
Expression Region	2-205
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.