



Recombinant Arabidopsis thaliana F-box protein At2g16450 (At2g16450)

Product Code	CSB-MP871080DOA
Abbreviation	At2g16450
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9SIV7
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MNPSPITIDL ILEILSRLPA KSVRRFHCVS KRWASIFGSP YFKELFLTRS STKPRLLFAI AEKGNKEKDC VWRFFSSPQL ENPYEKSSST LVATAEFHVR FSPDNLLICH YYDLKYFSIG YAFGLIYIG NRGRARPLIC NPTTGRYAIL PNRYTYRKAF SFFGFDPIDK QYKALSMVYP SGP GHSRVIT FGAGDLKWRR IKSLRHDIK SEGVCINGVL YYLGDTSDWS RVNGNHVTSG YMIVCFDVRS EKFTFIDVKR FCRLINYKKG LAVIYWEDDV DIQELYYKKG IDVEEYVENN VNADATNELC VWILADVEKQ EWSKHAYTWT DEKFFRRLVS IAGVTASGEI VFMSRKCNPQ QPFYVFYFNP ERNSLQRVEI QGFGEAVTKS CDVCTFVNHV EDLNVYDLKQ LKSVHPPLVE PEYYDSD
Source	Mammalian cell
Target Names	At2g16450
Protein Names	Recommended name: F-box protein At2g16450
Expression Region	1-427
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.