



Recombinant Arabidopsis thaliana Vacuolar protein sorting-associated protein 2 homolog 1 (VPS2.1)

| | |
|--------------------------|---|
| Product Code | CSB-EP871100DOA |
| Abbreviation | VPS2.1 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q9SKI2 |
| Product Type | Recombinant Protein |
| Immunogen Species | Arabidopsis thaliana (Mouse-ear cress) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MMNSIFGKRK TPAELLRENK RMLDKSIREI ERERQGLQTQ EKKLINEIKK TAKQQMGAV KVMAKDLIRT RHQIEKFYKL KSQLQGVSLR IQTLKSTQAM GEAMKGVTKA MGQMNRMNL PSLQKIMQEF ERQNEKMEMV SEVMGDAIDD ALEGDEEEEE TEDLVSQVLD EIGIDINQEL VNAPSGAVAV PAAKNKVVA EATGAEDSGG IDSDLQARLD NLRKM |
| Source | E.coli |
| Target Names | VPS2.1 |
| Protein Names | Recommended name: Vacuolar protein sorting-associated protein 2 homolog 1 Short name= AtVPS2-1 Alternative name(s): Charged multivesicular body protein 2 homolog 1 ESCRT-III complex subunit VPS2 homolog 1 SNF7-like protein |
| Expression Region | 1-225 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life |



of lyophilized form is 12 months at -20°C/-80°C.