



Recombinant Arabidopsis thaliana Ras-related protein RABF2b (RABF2B)

Product Code	CSB-MP871121DOA
Abbreviation	RABF2B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9SN68
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MAAAGNKSIN AKLVLLGDVG AGKSSLVLRV VKDQFVEFQE STIGAAFFSQ TLAVNDATVK FEIWDTAGQE RYHSLAPMYR RGAAAIIIVF DVTNQASFER AKKWWQELQA QGNPNMVMAL AGNKSDLLDA RKVTAEDAQT YQENGLFFM ETSAKTATNV KEIFYEIARR LPRVQPTENP TGMVLPDRAM DRAVSSSCCA
Source	Mammalian cell
Target Names	RABF2B
Protein Names	Recommended name: Ras-related protein RABF2b Short name= AtRABF2b Alternative name(s): Ras-related protein Ara-7 Ras-related protein Rab5B Short name= AtRab5B
Expression Region	1-200
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.