



Recombinant Arabidopsis thaliana Probable F-box protein At1g30780 (At1g30780)

Product Code	CSB-BP871203DOA
Abbreviation	At1g30780
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9SY17
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MASLSLVRLRY VVFLSFFILQ KPKETFSSRL QRSHAKVAQI GADDYPSSGR NTPVHDLGFP DFPSFQEAVA PPPPPDLPL LAPPLPDVPL LPPPAFPDFE KPRLPVPVWP SLPEYPPFPF VDQSAPYQGF APRFDESANW MPSTPSIPQV FPCLRFKSNL LRVEFWVSSG TAVKKGKDRFC NRPSKLDPNP VDLKTVRLPS KFRLKLLCR VESSEVCFGD GNEKDTNPSE IDLDSLPLFDL KMVILTRLISA KSLTNFKRVS KMWSSIIGSQ RFIDSFFTMS SKQSRCPKFP VNTRFVGYDP IDDQQKALSV SVPSRKRNL HKVLTGGGG QGWRHIEVTN APFSPVTGVV SIDGFVYYGA YSPTPPMNPV LVCFDVRSEK ISFIKAPNDV LQWGFDLWIL EDVERHEWSK QTCVFPSSVA WDVVVDIQMS FPGTNKAGIG DDEEFRRCSG FVDEGECHVR IAPQHVESIA RFKDPIMSRI LM
Source	Baculovirus
Target Names	At1g30780
Protein Names	Recommended name: Probable F-box protein At1g30780
Expression Region	1-482
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.