



Recombinant Arabidopsis thaliana Argininosuccinate synthase, chloroplastic (At4g24830)

Product Code	CSB-MP871222DOA
Abbreviation	At4g24830
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9SZX3
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	AVLSGDG TALTTDSKEA GLRGKLLKVV LAYSGGLDTS VIVPWLKENY GCEVVCFTAD VGQGIKELEG LEQKAKASGA SQLVVKDLTE EFVKDFIFPC LRAGAIYERK YLLGTSMARP VIAKAMVDVA AEVGADAVAH GCTGKGNDQV RFELTFFSLN PELKVVAPWR EWEIQGREDA IEYAKKHNP VPVTKKSIYS RDRNLWHLSH EGDLLDPAN EPKMDMYMMS VDPEDAPDQP EYIEIGIESG LPVALNGKAL SPATLLAELN TIGGKHGIGR IDMVENRLVG MKSRGVYETP GGTILFAAVQ ELESLLDRE SIQVKDTLAL KYAEMVYAGR WFDPLRESMD AFMEKITETT TGSVTLKLYK GSVSVTGRQS PNSLYRQDIS SFEGSEIYNQ ADAAGFIRLY GLPMKIRAML KKIS
Source	Mammalian cell
Target Names	At4g24830
Protein Names	Recommended name: Argininosuccinate synthase, chloroplastic EC= 6.3.4.5 Alternative name(s): Citrulline--aspartate ligase
Expression Region	74-494
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.