



Recombinant Dictyostelium discoideum Serine/threonine-protein phosphatase 6 catalytic subunit (ppp6c)

Product Code	CSB-EP871372DKK-B
Abbreviation	ppp6c
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9U9A3
Product Type	Recombinant Protein
Immunogen Species	Dictyostelium discoideum (Slime mold)
Purity	≥85% (SDS-PAGE)
Sequence	MSNALPLDEW VETARQCKYL PENDLKKLCE RVKELLLEES NVQPVRSPVT ICGDIHGQFY DLLELFKTGG EVPDTNYVFM GDFVDRGYYS LETFTYLLAL KARYPDKITL LRGNHESRQI TQVYGFYDEC QKQYGNVNAW KYCTSVDFDL TLAAIIDGKV LCVHGGGLSPK VRTLDAQIRII SRNLEIPHEG PFCDLMWSDP EDIEQWQPSP RGAGWLFSGK VTAEFEHING LNLICRAHQL VQEGYRYMFD NSLVTVWSAP NYCYRCGNVA SILSLNENLD RDFKIFQAVQ EERNIPTRPT MQYFF
Source	E.coli
Target Names	ppp6c
Protein Names	Recommended name: Serine/threonine-protein phosphatase 6 catalytic subunit Short name= PP6C EC= 3.1.3.16
Expression Region	1-305
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.