



Recombinant Human Zinc finger SWIM domain-containing protein 1 (ZSWIM1)

Product Code	CSB-MP871556HU
Abbreviation	ZSWIM1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BR11
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MLERLKAPWS AALQRKYFDL GIWTAPISPM ALTMLNGLLI KDSSPPMLLH QVNKTAQLDT FNYQSCFMQS VFDHFPEILF IHRTYNPRGK VLYTFLVDGP RVQLEGLHAR AVYFAIPAKE DTEGLAQMFQ VFKKFNPAWE RVCTILVDPH FLPLPILAME FPTAEVLLSA FHICKFLQAK FYQLSLERPV ERLLLTSLQS TMCSATAGNL RKLYTLLSNC IPPAKLPELH SHWLLNDRIW LAHRWRSRAE SSHYFQSLEV TTHILSQFFG TTPSEKQGMA SLFRYMQQNS ADKANFNQGL CAQNNHAPSD TIPESPKLEQ LVESHIQHSL NAICTGPAAQ LCLGELAVVQ KSTHLIGSGS EKMNIQILED THKVQPQPPA SCSCYFNQAF HLPSCRHILAM LSARRQVLQP DMLPAQWTAG CATSLDSILG SKWSETLDKH LAVTHLTEEV GQLLQHCTKE EFERRYSTLR ELADSWIGPY EQVQL
Source	Mammalian cell
Target Names	ZSWIM1
Protein Names	Recommended name: Zinc finger SWIM domain-containing protein 1
Expression Region	1-485
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.