



Recombinant Human SPRY domain-containing protein 5 (SPRYD5)

Product Code	CSB-EP871567HU-B
Abbreviation	SPRYD5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BSJ1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MNSGILQVFQ RALTCPICMN YFLDPVTIDC GHSFCRPCLY LNWQDTAVLA QCSECKKTR QRNLNTDIDL KNMAFIARKA SLRQFLSSEE QICGMHRETK KMFCEVDKSL LCLPCSNSQE HRNHIHCPIE WAAEERREEL LKKMQSLWEK ACENLRNLNM ETRTRCWKD YVSLRIEAIR AEYQKMPAFL HEEEQHHLER LRKEGEDIFQ QLNESKARME HSRELLRGM Y EDLKQMCHKA DVELLQAFGD ILHRYESLLL QVSEPVNPEL SAGPITGLLD SLSGFRVDFT LQPERANSHI FLCGDLRSMN VGCDPQDDPD ITGKSECFLV WGAQAFTSGK YYWEVHMGDS WNWAFGVCNN YWKEKRQNDK IDGEEGLFLL GCVKEDTHCS LFTTSPLVVQ YVPRPTSTVG LFLDCEGRTV SFVDVDQSSL YTIPNCSFS PPLRPIFCCS HF
Source	E.coli
Target Names	TRIM51
Protein Names	Recommended name: SPRY domain-containing protein 5 Alternative name(s): Tripartite motif-containing protein 51
Expression Region	1-452
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.