



Recombinant Human N-alpha-acetyltransferase 11 (NAA11)

Product Code	CSB-BP871569HU
Abbreviation	NAA11
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BSU3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MNIRNAQPDD LMNMQHCHNLL CLPENYQMKY YLYHGLSWPQ LSYIAEDEDG KIVGYVLAKM EEEPDDVPHG HITSLAVKRS HRRGLAQKL MDQASRAMIE NFNAKYVSLH VRKSNRPALH LYSNTLNFQI SEVEPKYYAD GEDAYAMKRD LSQMADELRR QMDLKKGGYV VLGSRENQET QGSTLSDSEE ACQQKNPATE ESGSDSKEPK ESVESTNVQD SSESSDSTS
Source	Baculovirus
Target Names	NAA11
Protein Names	Recommended name: N-alpha-acetyltransferase 11 EC= 2.3.1.88 Alternative name(s): N-terminal acetyltransferase complex ARD1 subunit homolog B Short name= hARD2 NatA catalytic subunit
Expression Region	1-229
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.