



Recombinant Human Leucine-rich repeat-containing protein 1 (LRRC1)

Product Code	CSB-BP871575HU
Abbreviation	LRRC1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BTT6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MFHCIPLWRC NRHVESIDKR HCSLVYVPEE IYRYARSLEE LLLDANQLRE LPEQFFQLVK LRKLGSLDNE IQRLPPEIAN FMQLVELDVS RNEIPEIPES ISFCKALQVA DFSGNPLTRL PESFPELQNL TCLSVNDISL QSLPENIGNL YNLASLELRE NLLTYLPDSL TQLRRLEELD LGNNEIYNLP ESIGALLHLK DLWLDGNQLS ELPQEIGNLK NLLCLDVSEN RLRLPPEEIS GLTSLTDLVI SQNLETIPD GIGKCLKLSI LKVDQNRLTQ LPEAVGECES LTELVLTENQ LLTLPKSIGK LKKLSNLNAD RNKLVSLPKE IGGCCSLTVF CVRDNRLTRI PAEVSQATEL HVLDVAGNRL LHLPLSLTAL KLKALWLSDN QSQPLLTFTQ DTDYTTGEKI LTCVLLPQLP SEPTCQENLP RCGALENLVN DVSDEAWNER AVNRVSAIRF VEDEKDEEDN ETRTLLRRAT PHPGELKHKM KTVENLRNDM NAAKGLDSNK NEVNHAIDRV TTSV
Source	Baculovirus
Target Names	LRRC1
Protein Names	Recommended name: Leucine-rich repeat-containing protein 1 Alternative name(s): LANO adapter protein LAP and no PDZ protein
Expression Region	1-524
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.