



Recombinant Human Multiple myeloma tumor-associated protein 2 (MMTAG2)

Product Code	CSB-EP871576HU
Abbreviation	MMTAG2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BU76
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MFGSSRGGVR GGQDQFNWED VKTDKQRENY LGNSLMAPVG RWQKGRDLTW YAKGRAPCAG PSREEELAAV REAEREALLA ALGYKNVKKQ PTGLSKEDFA EVCKREGGDP EEKGVDRLLG LGSASGSVGR VAMSREDKEA AKLGLSVFTH HRVESGGPGT SAASARRKPR AEDQTESSCE SHRKSKKEKK KKKKRKHKE KKKKDKHRR PAEATSSPTS PERPRHHHD SDSNSPCKR RKRGHSGDRR SPSRRWHDRG SEA
Source	E.coli
Target Names	MMTAG2
Protein Names	Recommended name: Multiple myeloma tumor-associated protein 2 Short name= hMMTAG2
Expression Region	1-263
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.