



Recombinant Human E3 ubiquitin-protein ligase TRIM62 (TRIM62)

Product Code	CSB-EP871583HU
Abbreviation	TRIM62
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BVG3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MACSLKDELL CSICLSIYQD PVSLGCEHYF CRRRCITEHWV RQEAQGARDC PECRRTFAEP ALAPSLKLAN IVERYSSFPL DAILNARRAA RPCQAHDKVK LFCLTDRALL CFFCDEPALH EQHQVTGIDD AFDELQRELK DQLQALQDSE REHTEALQLL KRQLAETKSS TKSLRTTIGE AFERLHRLLR ERQKAMLEEL EADTARTLTD IEQKVQRYSQ QLRKVGEGAQ ILQERLAETD RHTFLAGVAS LSERLKGKIH ETNLTYEDFP TSKYTGPLQY TIWKSFLQDI HPVPAALTLD PGTAHQRLIL SDDCTIVAYG NLHPQPLQDS PKRFDVEVSV LGSEAFSSGV HYWEVVVAEK TQWVIGLAHE AASRKGSIQI QPSRGFYCV MHDGNQYSAC TEPWTRLNVR DKLDKVG VFL DYDQGLLIFY NADDMSWLYT FREKFP GKLC SYFSPGQSHA NGKNVQPLRI NTVRI
Source	E.coli
Target Names	TRIM62
Protein Names	Recommended name: Tripartite motif-containing protein 62
Expression Region	1-475
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.