



Recombinant Human Complement C1q tumor necrosis factor-related protein 1 (C1QTNF1)

Product Code	CSB-MP871599HU
Abbreviation	C1QTNF1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BXJ1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	RVPHV QGEQQEWEGT EELPSPPDHA ERAEEQHEKY RPSQDQGLPA SRCLRCCDPG TSMYPATAVP QINITILKGE KGDRGDRGLQ GKYGKTGSAG ARGHTGPKGQ KGSMGAPGER CKSHYAAFSV GRKKPMHSNH YYQTVIFDTE FVNLYDHFNM FTGKFYCYVP GLYFFSLNVH TWNQKETYLH IMKNEEEVVI LFAQVGDRSI MQSQSLMLEL REQDQVWVRL YKGERENAIF SEELDTYITF SGYLVKHATE P
Source	Mammalian cell
Target Names	C1QTNF1
Protein Names	Recommended name: Complement C1q tumor necrosis factor-related protein 1 Alternative name(s): G protein-coupled receptor-interacting protein Short name= GIP
Expression Region	26-281
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.