



Recombinant Human Polycomb group RING finger protein 6 (PCGF6)

Product Code	CSB-BP871612HU
Abbreviation	PCGF6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BYE7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MEGVAVVTAG SVGAAKTEGA AALPPPPVVS PPALTPAPAA GEEGPAPLSE TGAPGCSGSR PPELEPERSL GRFRGRFEDE DEELEEEEEEL EEEEEEEEEED MSHFSLRLEG GRQDSEDEEE RLINLSELTP YILCSICKGY LIDATTITEC LHTFCKSCIV RHFYYSNRCP KCNIVVHQTQ PLYNIRLDRQ LQDIVYKLVI NLEEREKKQM HDFYKERGLE VPKPAVPQPV PSSKGRSKKV LESVFRIPPE LDMSLLLEFI GANEGTGHFK PLEKKFVRVS GEATIGHVEK FLRRKMGLDP ACQVDIICGD HLLQYQTLR EIRRAIGDAA MQDGLLVVHY GLVVSPLKIT
Source	Baculovirus
Target Names	PCGF6
Protein Names	Recommended name: Polycomb group RING finger protein 6 Alternative name(s): Mel18 and Bmi1-like RING finger RING finger protein 134
Expression Region	1-350
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This protein contains a RING finger motif, which is most closely related to those of polycomb group (PcG) proteins RNF110/MEL-18 and BMI1. PcG proteins are known to form protein complexes and function as transcription repressors. This protein has been shown to interact with some PcG proteins and act as a transcription repressor. The activity of this protein is found to be regulated by cell cycle dependent phosphorylation. Alternatively spliced transcript variants encoding different isoforms have been identified.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.