



# Recombinant Arabidopsis thaliana Probable pectate lyase 18 (At4g24780)

<b>Product Code</b>	CSB-EP871662DOA
<b>Abbreviation</b>	At4g24780
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9C5M8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	SPVPDP ESVVEEVHKS INASVAGRRK LGYLSCTTGN PIDDCWRCDP HWEQHRQRLA DCAIGFGKNA IGGRDGRIYV VTDSGNDNPV SPKPGTLRHA VVQDEPLWII FQRDMTIQLK EELIMNSFKT IDGRGASVHI SGGPCITIQY VTNIIHGIH IHDCKQGGNA MVRSSPRHFG WRTISDGDGV SIFGGSHVWV DHCSFSNCED GLIDAIMGST AITLSNNHMT HHDKVMLLGH SDTYSRDKNM QVTIAFNHFG EGLVQRMPRC RHGYFHVVNN DYTHWEMYAI GGSANPTINS QGNRFLAPNI RFSKEVTKHE DAPESEWKRW NWRSSGDLLE NGAFFTPSGG AASSYAKAS SLGAKPSSLV GPLTSTSGAL NCRKGSRC
<b>Source</b>	E.coli
<b>Target Names</b>	At4g24780
<b>Protein Names</b>	Recommended name: Probable pectate lyase 18 EC= 4.2.2.2 Alternative name(s): Pectate lyase A10
<b>Expression Region</b>	25-408
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.