



Recombinant Arabidopsis thaliana 14-3-3-like protein GF14 iota (GRF12)

Product Code	CSB-MP871665DOA
Abbreviation	GRF12
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9C5W6
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MSSSGSDKER ETFVYMAKLS EQAERYDEM V ETMKKVARVN SELTVEERNL LSVGYKNVIG ARRASWRIMS SIEQKEESKG NESNVKQIKG YRQKVEDELA NICQDILTII DQHLIPHATS GEATVFYK M KGDYRYLAE FKTEQERKEA AEQSLKGYEA ATQAASTE LP STHPIRLGLA LNFSVFY YEI MNSPERACHL AKQAFDEAIA ELDTLSEESY KDSTLIMQLL RDNLTLWTSD LPEDGGEDNI KTEESKQEQA KPADATEN
Source	Mammalian cell
Target Names	GRF12
Protein Names	Recommended name: 14-3-3-like protein GF14 iota Alternative name(s): General regulatory factor 12
Expression Region	1-268
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.