



# Recombinant Arabidopsis thaliana Pentatricopeptide repeat-containing protein At1g31790 (PCMP-A1)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-YP871673DOA   |
| <b>Abbreviation</b>      | PCMP-A1   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | Q9C6R9  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Arabidopsis thaliana (Mouse-ear cress)  |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | MKTLELALPP SPPSLVPSFN YNSTARSVGN DVRTNFDVQL FLRKPKHQKS<br>EPVVVIQQPQ IQPQNPSRC STSDILRLMD SLSLPGNEDI YSCLAKESAR<br>ENDQRGAHEL QVHIMKSSIR PTITFINRLL LMHVSCGRLD ITRQMFDRMP<br>HRDFHWAIV FLGCIEMGDY EDA AFLVSM LKHSQKGAFK IPSWILGCVL<br>KACAMIRDFE LGKQVHALCH KLG FIDEEDS YLSGSLIRFY GEFRLEDAN<br>LVLHQLSNAN TVAWAAKVTN DYREGEFQEV IRDFIEMGNH GIKKNVSVFS<br>NVLKACSWVS DGG RSGQQVH ANAIKLG FES DCLIRCLIE MYGKYGKVKD<br>AEKVKSSKD ETSVSCWNAM VASYMQNGIY IEAIKLLYQM KATGIKAHDT<br>LLNEAHLQM |
| <b>Source</b>            | Yeast   |
| <b>Target Names</b>      | PCMP-A1   |
| <b>Protein Names</b>     | Recommended name: Pentatricopeptide repeat-containing protein At1g31790   |
| <b>Expression Region</b> | 1-409   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.   |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life   |



of lyophilized form is 12 months at -20°C/-80°C.