



Recombinant Arabidopsis thaliana Serine carboxypeptidase-like 15 (SCPL15)

Product Code	CSB-BP871677DOA
Abbreviation	SCPL15
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9C7D2
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	SSIIRY LPGFEGPLPF ELETGYIGVG QKEEDQLFYY FIKSENNPEE DPLLWLTGG PGCSSFGLV YENGPLAFKV ETYNGSVPTL VSTTYSWTKV ANIIYLDQPV GTGFSYSRNP FADIPSDTGS VKRVNEFVRK WLA KHPEYFS NPFYVTGNSY SGKVIPAIQ EISNGNYICC KPQINLQGYV IGPNVAYYDH DKDFRIPFAH GVALISDEL FESLKASCGGS YSVVDPLNTE CLKLIEDYDK CVSGIYEELI LKSKCEHTSP DCYTYRYLLS EYWADNETVR RALKVVKGSK GTWERCDYRV LSNQDIKSSI PFHINNSIRG YRSLVISGDH DMTIPFLGTQ AWIRSLNYSI TEKWRPWMIL DQVAGYTKTY ANKMTLATVK GGGHTLEYKP EENSVLFKRW ISGQPL
Source	Baculovirus
Target Names	SCPL15
Protein Names	Recommended name: Serine carboxypeptidase-like 15 EC= 3.4.16.-
Expression Region	25-436
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.