



Recombinant Arabidopsis thaliana Putative F-box protein At1g33530 (At1g33530)

Product Code	CSB-BP871683DOA
Abbreviation	At1g33530
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9C800
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MTFSIEKLVD LKRKKWQSKL QNTKKQVLGR CLTNFSRLLS FRKKQETNGK SPRCSKATLA VVVHPTASTT HDDFLRKKQE RSGKSPSCSE TTLAVELPDV LVEEILQRLP VKYLVRLKSI SKGWKSLIES DHLAEKHLRL LEKKYGLKEI KITVERSTSK SICIKFFSRR SGMNAINSDS DLLLRVPGSC NGLVCVYELD SVYIYLLNPM TGVTRTLTPP RGTKLSVGFV IDVVTGTYKV MVLYGFDRVG TVVFDLDTNK WRQRYKTAGP MPLSCIPTPE RNPVFNGLSL FWLLASDFSE ILVMDLHTEK FRTLSPQNDM DDVDVSSGYI YMWSLEDRLC VSNVRQGLHS YVWVLVQDEL SEKWERTFRN LLGHVFPPLS LNSAWFSQTL VSPYQLSSST CIGSRQRQNS TSALFSRNGD TGAMDAPIEL HISVSTPMMP L
Source	Baculovirus
Target Names	At1g33530
Protein Names	Recommended name: Putative F-box protein At1g33530
Expression Region	1-441
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.