



Recombinant Mouse Mirror-image polydactyly gene 1 protein homolog (Mipol1)

Product Code	CSB-MP872021MO
Abbreviation	Mipol1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9D9F8
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MNSEQSIREL GNEVPSEDL E LPGRKTSNFQ VLQPCRDEGA SAECSIVECG KNYEPVISHQ VIPDLNKETS VAYLQKELEI LRASNTKLQE KLAKEDKEKR RLKLELQE KAAEADIAER TAALVEEVYF AQRERDEAIM CRLQLALKER DEAIAHVKHM EMSLKMLENI NPKENDMTLQ ELLDRINNAD TGIAIQKNGA VIVDTIYKTT GCRKRITAE E MSAVLEERDA ALSQCKQLHQ RLLDLKKQNN TSTNNTKHPT AKNNQEHTLK YNLNIASIN
Source	Mammalian cell
Target Names	Mipol1
Protein Names	Recommended name: Mirror-image polydactyly gene 1 protein homolog
Expression Region	1-279
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.