



Recombinant Mouse Calcium-binding protein 39-like (Cab39l)

Product Code	CSB-BP872039MO
Abbreviation	Cab39l
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9DB16
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MKKMPLFSKS HKNPAEIVKI LKDNLAILEK QDKKTDKASE EVSKSLQAMK EILCGTNDKE PPTEAVAQLA QELYSSGLLV TLIADLQLID FEGKKDVTQI FNNILRRQIG TRCPTVEYIS SHPHILFMLL KGYEAPQIAL RCGIMLRECI RHEPLAKIIL FSNQFRDFFK YVELSTFDIA SDAFATFKDL LTRHKVLVAD FLEQNYDTIF EDYEKLLQSE NYVTKRQSLK LLGELILDRH NFTIMTKYIS KPENLKLMMN LLRDKSPNIQ FEAHVFKVF VASPHKTQPI VEILLKNQPK LIEFLSSFQK ERTDDEQFAD EKNYLIKQIR DLKKAAP
Source	Baculovirus
Target Names	Cab39l
Protein Names	Recommended name: Calcium-binding protein 39-like Alternative name(s): MO25beta Mo25-like protein
Expression Region	1-337
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.