



Recombinant Mouse DnaJ homolog subfamily C member 10 (Dnajc10)

Product Code	CSB-YP872053MO
Abbreviation	Dnajc10
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9DC23
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	DQNFYSLL GVSKTASSRE IRQAFKKLAL KLHPDKNPNN PNAHGDFLKI NRAYEVLKDE DLRKKYDKYG EKGLEDNQGG QYESWSYYRY DFGIYDDDPE IITLERREFD AAVNSGELWF VNFYSPGCSH CHDLAPTWRE FAKEVDGLLR IGAVNCGDDR MLCRMKGVNS YPSLFIFRSG MAAVKYNGDR SKESLVAFAM QHVRSTVTEL STGNFVNAIE TAFAAGVGWL ITFCCKGEDC LTSQTRLRLS GMLDGLVNVG WVDCDAQDSL CKSLDTTAST TAYFPPGATL NDREKSSVLF LNSLDAKEIY MEIIHNLPDF ELLSANQLED RLAHHRWLVF FHFGKNENAN DPELKKLCTL LKNEHIQVGR FDCSSAPGIC SDLYVFQPCP AVFKGQGTKE YEIHHGKKIL YDILAFKES VNSHVTTLGP QNFPASDKEP WLVDFFAPWC PPCRALLPEL RCASTLLYGQ LKVGTLDCIT HEGLCNMYNI QAYPTTVVFN QSSIHEYEGH HSAEQILEFI EDLRNPSVVS LTPSTFNELV KQRKHDEVWM VDFYSPWCHP CQVLMPEWKR MARTLTGLIN VGSVDCQQYH SFCTQENVQR YPEIRFYPQK SSKAYQYHSY NGWNRDAYSL RSWGLGFLPQ ASIDLTPQTF NEKVLQGKTH WVVDYFAPWC GPCQNFAPFP ELLARMIKKG VRAGKVDCCA YPQTCQKAGI KAYPSVKLYQ YERAKKSIWE EQINSRDAKT IAALIYGKLE TLQSQVKNK DEL
Source	Yeast
Target Names	Dnajc10
Protein Names	Recommended name: DnaJ homolog subfamily C member 10 Alternative name(s): ER-resident protein ERdj5 Endoplasmic reticulum DnaJ-PDI fusion protein 1 J domain-containing protein disulfide isomerase-like protein Short name= J domain-co
Expression Region	33-793
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein



Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.