



Recombinant Mouse Guanine nucleotide-binding protein G (k) subunit alpha (Gnai3)

Product Code	CSB-BP872054MO
Abbreviation	Gnai3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9DC51
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	GCTLSAEDK AAVERSKMID RNLREDGEKA AKEVKLLLLL AGESGKSTIV KQMKIIHEDG YSEDECKQYK VVVYSNTIQS IIAIIRAMGR LKIDFGESAR ADDARQLFVL AGSAEEGVM T SELAGVIKRL WRDGGVQACF SRSREYQLND SASYLNDLD RISQNTNYIPT QQDVL RTRVK TTGIVETHFT FKELYFKMFD VGGQRSERKK WIHCFEGVTA IIFCV ALSDY DLVLAEDEEM NRMHESMKLF DSICNNKWFT DTSIILFLNK KDLFEEKIKR SPLTICYPEY TGSNTYEEAA AYIQCQFEDL NRRKDTKEVY THFTCATDTK NVQFV FDAVT DVIKNNLKE CGLY
Source	Baculovirus
Target Names	Gnai3
Protein Names	Recommended name: Guanine nucleotide-binding protein G(k) subunit alpha Alternative name(s): G(i) alpha-3
Expression Region	2-354
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.