



Recombinant Human Guanine nucleotide-binding protein-like 1 (GNL1)

Product Code	CSB-EP009623HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P36915
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	<p>MPRKKPFSVK QKKKQLQDKR ERKRGLQDGL RSSSNSRSGS RERREEQTD TSDGESVTHHI RRLNQQPSQG LGPRGYDPNR YRLHFERDSR EEVERRKRAA REQVLQPVSA ELLELDIREV YQPGSVLDFP RRPPWSYEMS KEQLMSQEER SFQDYLKGIH GAYSSEKLSY FEHNLETWRQ LWRVLEMSDI VLLITDIRHP VVNFPPALYE YVTGELGLAL VLVLNKVDLA PPALVVAWKH YFHQHYPLH VVLFSTFPRD PRTPQDPSSV LKKSRRRGRG WTRALGPEQL LRACEAITVG KVDLSSWREK IARDVAGATW GNGSGEEEEEE EDGPAVLVEQ QTDSAMEPTG PTQERYKDG VVTIGCVGFNP VGKSSLINGL VGRKVVSVSR TPGHTRYFQT YFLTPSVKLC DCPGLIFPSL LPRQLQVLG IYPIAQIQEP YTAVGYLASR IPVQALLHLR HPEAEDPSAE HPWCAWDICE AWA EKRGYKT AKAARNDVYR AANSLRLAV DGRSLCFHP PGYSEQKGTW ESHPETTEL V VLQGRVGPAG DEEEEEEEEL SSSCEEEGEE DRDADEEGEG DEETPTSAPG SSLAGRNPYA LLGEDEC</p>
Source	E.coli
Target Names	GNL1
Protein Names	Recommended name: Guanine nucleotide-binding protein-like 1 Alternative name(s): GTP-binding protein HSR1
Expression Region	1-607
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Target Details	The GNL1 gene, identified in the human major histocompatibility complex class I region, shows a high degree of similarity with its mouse counterpart. The GNL1 gene is located less than 2 kb centromeric to HLA-E, in the same transcriptional orientation. GNL1 is telomeric to HLA-B and HLA-C.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.