



Recombinant Tetraodon fluviatilis Hsp90 co-chaperone Cdc37 (cdc37)

Product Code	CSB-BP872081TLV
Abbreviation	cdc37
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9DGQ7
Product Type	Recombinant Protein
Immunogen Species	Tetraodon fluviatilis (Green pufferfish) (Chelonodon fluviatilis)
Purity	>85% (SDS-PAGE)
Sequence	MSRIDYSVWD HIEVSDDDEDV SHPNIDTPSL FRWRHQARVE RMEDFKKKGD DLNKGLQECR RKLAEAQKKV RELSISAAGD AKAELTKTQT EEKCLKKEER DWARDLEEHN HEEKKMPWNV DTLKDGFSK SIVNVKADSA EDTEEEKEKK HKTFVERYEK QIKHFGMLHR WDDSQRYLSD NPDLVCEETA NYLVIMCIDL EVEEKHALME QVAHQTIVMQ FILELAKSLK VDPRGCFRQF FEKIKTADQQ YQDAFNDELE SFKERVRGRA KIRIEKAMKE YEEEEERQKRL GPGGLDPVEV YESLPPMQK CFDEKDIQML QDVITKMDPT EAKAHMKRCI DSGLWVPNSK MDESDEKEDG STYEEVKQEQ EEPAPKE
Source	Baculovirus
Target Names	cdc37
Protein Names	Recommended name: Hsp90 co-chaperone Cdc37 Alternative name(s): Hsp90 chaperone protein kinase-targeting subunit p50Cdc37
Expression Region	1-377
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.