



# Recombinant Rat Alkyldihydroxyacetonephosphate synthase, peroxisomal (Agps)

<b>Product Code</b>	CSB-MP872120RA
<b>Abbreviation</b>	Agps
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9EQR2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	KARRAA SAAGASPAAS PAAPESGTIP KKRQELMKWN GWGYNDSKFL LNKKGQVELT GKRYPLSGLA LPTLKDWIQN TLGVNLEHKT TSKPSINPSE APPSIVNEDF LQELKEAHIS YSQDADDRVF RAHGHCLHEI FLLREGMFER IPDIVVWPTC HDDVVKIVNL ACKYNLCIIP IGGGTSVSYG LMCPAETRTR IISLDTSQMN RILWVDENNL TAHVEAGITG QDLERQLKES GYCTGHEPDS LEFSIVGGWI STRASGMKKN VYGNIEDLVV HMKMVTPRGV IEKSSQGPRM STGPDIIHFI MGSEGLGVI TEATIKIRPT PEYQKYGSVA FPNFEQGVAC LREIAKQRCA PASIRLMDNQ QFQFGHALKP QVSSIFTSFL DGFKKFYITK FKGFDPNQIS VATLLFEGDR EKVLQHEKQV YDIAAKFGGL AAGEDNGQRG YLLTYVIAYM RDLGLEYYVV GESFETSAPW DRVIDLCRNV KERIRRECKE RGVQFAPLST CRVTQTYDAG ACIYFYFAFN YRGISDPLTV FEQTEAAARD EILANGGSL S HHHGVGKLRK QWLKESISDV GFGMLKSVKD YVDPSNIFGN RNLL
<b>Source</b>	Mammalian cell
<b>Target Names</b>	Agps
<b>Protein Names</b>	Recommended name: Alkyldihydroxyacetonephosphate synthase, peroxisomal Short name= Alkyl-DHAP synthase EC= 2.5.1.26 Alternative name(s): Alkylglycerone-phosphate synthase
<b>Expression Region</b>	45-644
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This gene is a member of the FAD-binding oxidoreductase/transferase type 4 family. It encodes a protein that catalyzes the second step of ether lipid biosynthesis in which acyl-dihydroxyacetonephosphate (DHAP) is converted to alkyl-DHAP by the addition of a long chain alcohol and the removal of a long-



chain acid anion. The protein is localized to the inner aspect of the peroxisomal membrane and requires FAD as a cofactor. Mutations in this gene have been associated with rhizomelic chondrodysplasia punctata, type 3 and Zellweger syndrome.

---

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.