



# Recombinant Mouse Torsin-3A (Tor3a)

<b>Product Code</b>	CSB-MP872125MO
<b>Abbreviation</b>	Tor3a
<b>Storage</b>	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
<b>Uniprot No.</b>	Q9ER38
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	QEADPTPW PSVKGLKEQL RKAGALSKRY WELFCTLWP DHCEDQETPV PPLGWSLPLW GRRSLDVLTA WLCRFQDCCS GGGDCRISNN LTGLESDLRV RLHGQHLASK LVLRAVKGYL EMPQVGKALA LSFHGWSGTG KNFLARILMD NLYRDGMRS D CVKMFISTFH FPHPKYVDY KEELQRQMQE TQWRCHQSTF VFDEAEKLHP GLELLEPYL EPRSPEARV EAPRAIFLFL SNLGGSVINE VVLSELLKAGW SREEITTQHL EVPLQAEIME AADSSFGSSG LLKKHLIDHF IPFLPLEYCH VRLCVRDAFL GQDLPYTEET LDEIAKMMTY VP EEERLFSS QGCKSISQRI NLFLP
<b>Source</b>	Mammalian cell
<b>Target Names</b>	Tor3a
<b>Protein Names</b>	Recommended name: Torsin-3A Alternative name(s): ATP-dependent interferon-responsive protein Torsin family 3 member A
<b>Expression Region</b>	22-385
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>