



Recombinant Arabidopsis thaliana F-box protein SKIP1 (SKIP1)

Product Code	CSB-MP872199DOA
Abbreviation	SKIP1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FDX1
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MEEDGSDWGG LAPEILINII SRLTIQELWT GPMFVQKSWL TVCRDPYLWS IFDLEPWFDY YPESTHLWSP EFEQKVDLML RSVVDWSEGG LTKIRVRHCS DHALSYAADR CPNLQVLAIR SSPNVTDASM TKIAFRCSRSL KELDYSYCHE ISHDTLVMIG RNCPNLRILK RNLMDWSSRH IGSVPTEYLD ACPQDGDTEA DAIGKHMNL EHLEIQFSRL SVKGLASICE GCPKLEYLDL FGCVHLSSRD ITSNVSRLKW LKEVKKPDVY VPRSGDVAQT ERYGHWRLYD ERFDIQAMRI
Source	Mammalian cell
Target Names	SKIP1
Protein Names	Recommended name: F-box protein SKIP1 Alternative name(s): SKP1-interacting partner 1
Expression Region	1-300
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.