



Recombinant Arabidopsis thaliana Fasciclin-like arabinogalactan protein 13 (FLA13)

Product Code	CSB-YP872209DOA
Abbreviation	FLA13
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FFH6
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	APAPG PAGPINITAI LEKGGQFVTL IRLNNTTQIG NQINIQINSS SEGMTVLAPT DNAFQNLKPG TLNKLSPDDQ VKLILYHVSP KFYTLEDLLS VSNPVRTQAS GRDVGGVYGL NFTGQGNQVN VSTGVVETRL STSLRQERPL AVYVVDMVLL PEEMFGERKI SPMAPPKSK SPDVSDDSES SKKAAAPSES EKSG
Source	Yeast
Target Names	FLA13
Protein Names	Recommended name: Fasciclin-like arabinogalactan protein 13
Expression Region	26-224
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.