



Recombinant Arabidopsis thaliana Probable F-box protein At5g36000 (At5g36000)

Product Code	CSB-EP872217DOA-B
Abbreviation	At5g36000
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FGB8
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MNTRSGDAEG DIRGKMIAPV RDGNGGQKRK LVQSNDIQRD EDGGAKRRII QSSDQKNGKI LRGIHGCVSP RCSAPTYQSR FSWYEQDIWT YITRFLDGKS LVKLGATNKW FYKIAMEDTV WRFACLRLDQ VPETFPVSST WIKIYASAFD GSHSYLFHQK EKHIDWMRLG AFVLDSRTSF LTESLSGRLK VPTEGTIEQM LQSSGSLIN DIKSGIWIAD LQLLRCPLCD LSTCDGTMQT LDVRHIELFL SEGYKDGSDW YNLIGSHKLE KDASAACGAI FDLKHLKSS SSGILNLKSW TGEADDSQPK AVIAPHAVAV HTRLQQNEGI LVKYHTMKAG TDGDIVSIRI SQLL
Source	E.coli
Target Names	At5g36000
Protein Names	Recommended name: Probable F-box protein At5g36000
Expression Region	1-355
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.