



Recombinant Arabidopsis thaliana Putative F-box protein At5g42430 (At5g42430)

| | |
|--------------------------|---|
| Product Code | CSB-MP872241DOA |
| Abbreviation | At5g42430 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q9FIH4 |
| Product Type | Recombinant Protein |
| Immunogen Species | Arabidopsis thaliana (Mouse-ear cress) |
| Purity | ≥85% (SDS-PAGE) |
| Sequence | MNGEENS DSI PIDLILEILS RLPKASITRF HCVSKLWGSM LCRPYFNELF LTISSAR PRL LFAFSKHGEW RFFSSPQPQN PYGKSSFVAT ADFHTKFSQN LNICNYT SGL VYFSAMWITK ADVICNPSTG HYAMPLKLLL TYGETRSFFL FDPVGKQFKV LLMNKINNNE TKDIHILTLG TRKVRWRKIQ QCPLIHIVSH EWICINGALY YIAYNIDDFL GYIVCFD VRS EKFKCLNLNQ DCFSERSTKL IYYKGLGVV NLKYAHGGGF PLKLCMWVLE DVEKQEWTTTS VYTLRDEDRV VKVYYDLFIV GITATGEIVL AKKKVCKPFY VFYFNLERNT LLSVEIQGFG EYQSCCSVHA FVDHVEDLNV YAFVEHMKKT YEYDATSISP SEQKL |
| Source | Mammalian cell |
| Target Names | At5g42430 |
| Protein Names | Recommended name: Putative F-box protein At5g42430 |
| Expression Region | 1-395 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |