



Recombinant Arabidopsis thaliana Quercetin 3-O-methyltransferase 1 (OMT1)

Product Code	CSB-EP872258DOA
Abbreviation	OMT1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FK25
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MGSTAETQLT PVQVTDDEAA LFAMQLASAS VLPMALKSAL ELDLLEIMAK NGSPMSPTTEI ASKLPTKNPE APVMLDRILR LLTSYSVLTC SNRKLSGDGV ERIYGLGPVC KYLTKNEDGV SIAALCLMNQ DKVLMESWYH LKDAILDGGI PFNKAYGMSA FEYHGTDPHF NKVFNNGMSN HSTITMKKIL ETYKGFEGLT SLVDVGGGIG ATLMIVSKY PNLKGINFDL PHVIEDAPSH PGIEHVGGDM FVSVPKGDAL FMKWICHDWS DEHCVKFLKN CYESLPEDGK VILAECILPE TPDSSLSTKQ VVHVDCIMLA HNPGGKERTE KEFEALAKAS GFKGIKVVCD AFGVNLIELL KKL
Source	E.coli
Target Names	OMT1
Protein Names	Recommended name: Quercetin 3-O-methyltransferase 1 Short name= AtOMT1 EC= 2.1.1.76 Alternative name(s): Flavonol 3-O-methyltransferase 1
Expression Region	1-363
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.