



Recombinant Arabidopsis thaliana Probable MYST-like histone acetyltransferase 1 (HAG4)

Product Code	CSB-BP872276DOA
Abbreviation	HAG4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FLF7
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MGSSADTETA MIIATPASNH NNPATNGGDA NQNHTSGAIL ALTNSESDAS KKRRMGVLPL EVGTRVMCQW RDGKYHPVKV IERRKNYNGG HNDYEYYVHY TEFNRRLDEW IKLEQLDLDS VECALDEKVE DKVTSCLKMTR HQRKIDETH VEGHEELDAA SLREHEEFTK VKNIATIELG KYEIETWYFS PFPPEYNDCV KLFFCEFCLS FMKRKEQLQR HMRKCDLKHG PGDEIYRSST LSMFEVDGKK NKVYAQNLCY LAKLFLDHTK LYDVLDFLF YILCECDRGG CHMVGYSKE KHSEEAYNLA CILTLPPYQR KGYGKFLIAF SYELSKKEGK VGTPERPLSD LGLVSYRGYW TRILLDILKK HKNISIKEL SDMTAIKAED ILSTLQSELE IQYRKGQHVI CADPKVLDRH LKAAGRGLD VDVSKMIWTP YKEQS
Source	Baculovirus
Target Names	HAM1
Protein Names	Recommended name: Probable MYST-like histone acetyltransferase 1 EC=2.3.1.48
Expression Region	1-445
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.