



Recombinant Arabidopsis thaliana Putative F-box/kelch-repeat protein At5g24040 (At5g24040)

Product Code	CSB-EP872282DOA
Abbreviation	At5g24040
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FLV8
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MVKWSELPPE ILHLISLKID NPFDLIHFRS VCSFWRSSSL LKFRHMTSLR CPLPLDPGGC GDDCHILSSR VYLLKCPNRD RPQYWMFKLQ AKENGEVVLH SLFLRRNSSA YGCLYPSLSI DLLNCQIFEL AQEHVACYSE WFELFECISK CEERIGFMGL NREKNEYMIL GKLSFNGLAM YRSVDNRWTE LEITPDSFFE GIVPFKGFY AIDRTGKTTV VDPTLEVNTF QRSRPCDKTR KRWLLMTGDK LILVEMCTS RYDFHIPNIR EKKIWFSE LNEERNDWDQ VEDMDGRVLF LEHYCTFSCL ATEIPGFRAN SIIFMDLWGG SNSYELESIL VYEFNEKGVR SLRDKLEYIE LFPFPPGWVI SNG
Source	E.coli
Target Names	At5g24040
Protein Names	Recommended name: Putative F-box/kelch-repeat protein At5g24040
Expression Region	1-373
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.