



Recombinant Solanum tuberosum Isovaleryl-CoA dehydrogenase 2, mitochondrial (IVD2)

Product Code	CSB-BP872315FIG
Abbreviation	IVD2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FS87
Product Type	Recombinant Protein
Immunogen Species	Solanum tuberosum (Potato)
Purity	≥85% (SDS-PAGE)
Sequence	ST SLLFDDTQKQ FKESVAQFAQ ENIAPHAEKI DRTNYFPQDV NLWKLGMGNFN LLGITVPEEY GGLGLGYLYH CIAMEEISRA SGSVGLSYGA HTNLCINQLV RNGTHEQKQK YLPKLISGEH VGALAMSEPN AGSDVVS MKC KADRVEGGYV LNGNKMWCTN GPTAQTLLVY AKTDVTAGSK GITAFIIEKG MTGFSTAQKL DKLGMRGSDT CELVFENCFV PEENVLGQVG RGVYVLSGL DLERLVLASG PVGIMQACLD VVLPYVKQRE QFGRPIGEFQ FVQGKVADMY TSMQSSRSYL YSVARECD SG TINTKDCAGV ILSAAERATQ VALQAIQCLG GNGYVNEYPT GRFLRDAKLY EIGAGTSEIR RMIIGRELFK EQ
Source	Baculovirus
Target Names	IVD
Protein Names	Recommended name: Isovaleryl-CoA dehydrogenase 2, mitochondrial Short name= IVD 2 EC= 1.3.99.10
Expression Region	29-412
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.